



Validation Report for SARS-CoV-2 PriSt MTM

CA-VALRPT-LAB-003
Version number 2.2

CA-VALRPT-LAB-003 Validation Report for SARS-CoV-2 PriSt MTM

Copy of version 2.2 (approved and current)

Last Approval or
Periodic Review Completed 2/19/2021

Next Periodic Review
Needed On or Before 8/19/2021

Effective Date 3/2/2021

Uncontrolled Copy printed on 5/14/2021 12:06
PM

Printed By Paven Gahir

Organization PerkinElmer Genomics California

Author
Lynn Deng

Comments for version 2.0 (last major revision)

Added Section 2.3

Comments for version 2.2 (this revision)

Revision history corrected.

Approval and Periodic Review Signatures

Type	Description	Date	Version	Performed By	Notes
Approval	Lab Director	2/19/2021	2.0	Adam Rosendorff	
Approval	Approval by QA Director	2/19/2021	2.0	<i>Lora J. H. Bean</i> Lora Bean	
Approval	Lab Director	11/23/2020	1.0	<i>Shantelle Lucas</i> Shantelle Lucas	

Version History

Version	Status	Type	Date Added	Date Effective	Date Retired
2.2	Approved and Current	Minor revision	3/2/2021	3/2/2021	Indefinite
2.2	Retired	Minor revision	3/2/2021	3/2/2021	3/2/2021
2.0	Retired	Major revision	2/19/2021	2/19/2021	3/2/2021
1.0	Retired	Initial version	11/21/2020	11/23/2020	2/19/2021

Approved: 2/19/2021 by Adam Rosendorff



Validation Summary for SARS-CoV-2 Testing using PrimeStoreMTM

CA-VALRPT-CYT-2001

1. PURPOSE

- 1.1 To provide additional data to validate a molecular diagnostic assay that utilizes real-time polymerase chain reaction (PCR) in order to determine the presences or absences of SARS-CoV-2 RNA in human nasal swabs at CDPH Branch Laboratory, Valencia, CA.

2. SCOPE

- 2.1 This procedure applies to individuals performing the SARS-CoV-2 Assay at CDPH Branch Laboratory, Valencia, CA.
- 2.2 This validation report is for SARS-CoV-2 RT-qPCR Reagent kit (Cat # 2019-nCoV-PCR-AUS) and is a supplement to CA-VALRPT-LAB-002.
- 2.3 This SARS-CoV-2 Assay is approved for each non-waived test for clinical use

3. DEFINITIONS

- 3.1 IC: An exogenous internal control comprised of TE buffer and bacteriophage MS2, provided in the kit (Cat #: 2019-nCoV-PCR-AUS, tube: nCoV Positive Control)
- 3.2 Ct or Cq: Cycle number at which PCR response starts to become exponential
- 3.3 AN: Anterior Nasal
- 3.4 NSP: Nasopharyngeal
- 3.5 ORF1ab and N genes: 2 different coding areas in the SARS-CoV2 genome interrogated in this validation using RT-PCR
- 3.6 FAM, ROX, HEX/VIC: Abbreviations for 3 fluorescent dyes used in reporter oligonucleotides in RT-PCR assay. FAM probe reports on SARS-Cov2 N gene, ROX on ORF1ab gene, and HEX reports on the IC (bacteriophage MS2)
- 3.7 RT-PCR: real-time PCR
- 3.8 SARS-CoV-2: Severe acute respiratory syndrome-related coronavirus of the genus Betacoronavirus, strain 2
- 3.9 PrimerStoreMTM – inactivating molecular transport media used in these studies.

4. ROLES AND RESPONSIBILITIES

- 4.1 It is the responsibility of the laboratory and delegates to establish any methods required to perform this assay. This must be properly documented in a validation plan, report, standard operating procedures, and any other required documentation.



Validation Summary for SARS-CoV-2 Testing using PrimeStoreMTM

CA-VALRPT-CYT-2001

004)

- 6.1.6 SARS-CoV-2 RT-PCR Set-up Using Janus G3 (CA-PCR-SOP-001)
- 6.1.7 SARS-CoV-2 RT-PCR Using the Analytik Jena (CA-PCR-SOP-002)
- 6.1.8 SARS-CoV-2 Assay Data Extraction (CA-RPT-SOP-001)
- 6.1.9 Analysis and Reporting of SARS-CoV-2 Assay (CA-RPT-SOP-002)

6.2 Equipment and Systems with ID

Table 1. Equipment and System IDs					
Equipment	Manufacturer	Model	Part/Cat#	Purpose	Environmental Requirements Temperature/ Humidity
chemagic™	PerkinElmer®	360 or 360-D	20240056	DNA/RNA Extraction	Temp: 18 - 35 °C Humidity: < 80 %
Biosafety Hood	Any	Any	Any	Prepare master mix	Dependent on model
JANUS G3	PerkinElmer®	Any	Any	Sample Reformatting and PCR workstation	Temp: 15 – 35 °C Humidity: 60 – 80 %
Microcentrifuge	Any	Any	Any	Collect samples at bottom of well	Dependent on model
Vortex mixer	Any	Any	Any	Mix reagents and/or samples	Dependent on model
Plate Centrifuge	Any	Any	Any	Collect samples at bottom of well	Dependent on model
Pipettes (single and multi-channel) p10, p200, p1000	Any	Any	Any	Pipette reagent and samples	Dependent on model



Validation Summary for SARS-CoV-2 Testing using PrimeStoreMTM

CA-VALRPT-CYT-2001

Table 1. Equipment and System IDs					
Equipment	Manufacturer	Model	Part/Cat#	Purpose	Environmental Requirements Temperature/ Humidity
Heratherm Oven	Thermo Scientific	Any	Any	Heat samples to kill virus	Temp: 15 - 31 °C Humidity: 10 - 85 %
qPCR Real-Time System	Analytica Jena	AJ384	844-00569-4	RT-PCR	Temp: 15 - 31 °C Humidity: 10 - 85 %

6.3 Reagents, Supplies, and Materials

Table 2. Reagents					
Reagent	Vendor	Manufacturer	Part/Cat#	Purpose	Storage Conditions
chemagic™ Viral DNA/RNA 300 Kit H96 component <u>Magnetic Beads B</u> 1 bottle, 150 mL Prep: Ready to use Exp: Unopened, according to the labeled Opened: 30 days	PerkinElmer®	AUS	CMG-1033-S	Bind RNA/DNA	2 - 25 °C
chemagic™ Viral DNA/RNA 300 Kit H96 component <u>Lysis Buffer 1</u> 1 bottle, 320 mL Prep: Ready to use Exp: Unopened, according to the label Opened: 30 days	PerkinElmer®	AUS	CMG-1033-S	Lyses cells or other DNA source to get the RNA/DNA into solution	2 - 25 °C



Validation Summary for SARS-CoV-2 Testing using PrimeStoreMTM

CA-VALRPT-CYT-2001

Table 2. Reagents					
Reagent	Vendor	Manufacturer	Part/ Cat#	Purpose	Storage Conditions
chemagic™ Viral DNA/RNA 300 Kit H96 component Poly (A) RNA Buffer Mix 10 bottles, 2 mL Exp: Unopened, according to label Opened: 30 days.	PerkinElmer®	AUS	CMG-1033-S	Added to Poly (A) RNA to activate.	15 - 25 °C
PrimeStoreMTM – Inactivating molecular transport media	Any	Longhorn vaccines & diagnostics LLC	Variable	Media to create contrived LOD samples	15 - 25 °C
AccuPlex™ SARS-CoV-2 Molecular Controls Kit – Full Genome	SeraCare	SeraCare Life Science, Inc.	0505-0159	Used to create contrived samples	2 - 8 °C
New Coronavirus Nucleic Acid Detection Kit	PerkinElmer®	PerkinElmer®	2019-nCoV-PCR-AUS	RT-PCR kit	-25 to -15°C

Table 3. Materials and Supplies					
Material and Supplies	Vendor	Manufacturer	Part/Cat #	Purpose	Storage Conditions
2 mL Deep Well Plate	PerkinElmer®	PerkinElmer®	CMG-555	Load sample/lysis and Elution Buffer B onto the chemagic™ 360	15 - 25 °C
96 Rod Head Disposable Tips	PerkinElmer®	PerkinElmer®	49017-0006	Protect 96 Rod Magnetic Head (on chemagic™ 360) from contamination	15 - 25 °C

Validation Summary for SARS-CoV-2 Testing using PrimeStoreMTM

CA-VALRPT-CYT-2001

Table 3. Materials and Supplies					
Material and Supplies	Vendor	Manufacturer	Part/Cat #	Purpose	Storage Conditions
Low-well Plate	PerkinElmer®	PerkinElmer®	CMG-555-1	Load Magnetic Beads B onto the chemagic™ 360	15 - 25 °C
Heat Sealing Foil	Approved vendor	Approved vendor	Approved vendor	Seal plates	N/A
Reagent trough	Any	Any	Any	Preserve reagent stock bottle and allow for multi-channel pipetting	N/A
Tips (10 µL, 200 µL, 1000 µL) compatible with pipettes	Any	Any	Any	Pipette reagent and samples	Dependent on brand
384 Plates for RT-PCR	Approved vendor	Approved vendor	Approved vendor	Hold samples and master mix	N/A
qPCR Seal	Approved vendor	Approved vendor	Approved vendor	Seal 384 well plate	N/A; Expiration Date is written on the box
1.7 mL to 50 mL tube	Any	Any	Any	For master mix preparation	Dependent on brand

6.4 Training Requirements

6.4.1 Personnel that performing this assay will require proper training. This training will include complete review and understanding of all associated documentation and hands-on laboratory training in the performance of this assay. This training will be documented on the Process Training Form.



Validation Summary for SARS-CoV-2 Testing using PrimeStoreMTM

CA-VALRPT-CYT-2001

7. TESTING PERFORMED

7.1 Limit of Detection (LOD): LOD Range Finding

7.1.1 Contrived samples were created in triplicate. The LODs tested were 10, 20, 40, 60, 120, 180 copies per mL in a total volume of 1 mL per sample. In addition, samples with 0 copy were tested.

7.1.2 Samples were tested both with and without heat inactivation.

7.1.3 A total of 36 samples were run:

7.1.3.1 21 without heat: each LOD in triplicate

7.1.3.2 23 with heat: each LOD in triplicate except 0 copy in 5 replicates.

7.2 Confirmation of LOD

7.2.1 The LOD Confirmation run established the LOD at 120 copies / mL. Samples were created in 20 replicates at 120 copies mL as 1X, and in a level below and above 60 copies / mL and 180 copies /mL.

7.2.2 Due to limited availability of clinical samples, LOD confirmation was performed in PrimeStoreMTM rather than spiked previously tested negative samples.

7.3 Clinical Evaluation

7.3.1 60 previously tested clinical samples (30 samples previously tested as negative and 30 previously tested as positive) were tested. These clinical samples were obtained from a San Mateo County Public Health (CLIA 05D0857622) Laboratory performing SARS-CoV-2 testing.

7.3.2 All previously tested samples were collected in VTM and combined with PrimeStoreMTM in a 1:1 ratio for this study.

7.4 Verification of precision and accuracy

7.4.1 Inter-run Precision: 60 samples were run on different equipments by different technologists.

7.4.2 Intra-run Precision: 12 samples were run in duplicate.

7.4.3 Accuracy: Determined by comparison to previously tested clinical samples.

8. RESULTS

8.1 LOD Finding (See attached Excel sheet)



Validation Summary for SARS-CoV-2 Testing using PrimeStoreMTM

CA-VALRPT-CYT-2001

11. APPENDICES

11.1 Appendix A: Evidence of Dr. Lucas’s Approval

12. REVISION HISTORY

Version	Summary of Changes	Date
2.2	Appendix A Added & Revision History updated	March 2, 2021
2.1	Inadvertently removed – no changes lost	
2	Section 2.3 Added	Feb 19, 2021
1	New Document	

APPENDIX A

Evidence of Dr. Lucas’s approval on 24Oct2020 t 1:55pm is evidenced by her email to LFS

From: [Shantelle Lucas](#)
To: [Bean, Lora](#)
Cc: [Dowless Kozar, Holli](#)
Subject: [External] FW: CDPH Branch Laboratory Validation
Date: Wednesday, February 24, 2021 11:58:43 AM
Attachments: [PrimeStoreMTM_10_24_2020.zip](#)

Use caution when opening links or attachments.

Lora,

See attached.

Shantelle

From: Shantelle Lucas
Sent: Saturday, October 24, 2020 1:55 PM
To: Eleco, Elsa@CDPH <Elsa.Eleco@cdph.ca.gov>; Flores, Elaine@CDPH <Elaine.Flores@cdph.ca.gov>; Thomas, Robert@CDPH <Robert.Thomas@cdph.ca.gov>; Haleh Farzanmehr <hfarzanmehr@genexlpc.com>; Farzanmehr, Haleh@CDPH <Haleh.Farzanmehr@cdph.ca.gov>
Subject: CDPH Branch Laboratory Validation

Hello Elsa, Elaine and Bob,

Attached are the documents for the CDPH Branch Laboratory Validation . Of note, samples only collected in MTM will be used at this time for testing. If you have any questions, please let us know.

Thank you,
Shantelle Lucas

Uncontrolled Copy