

## Micro Learning Video: CCC Changing an Adequacy Code

In this video, you will learn how to change adequacy codes.

1. Click the Accession # of the specimen record.
2. Click the Pregnancy Information tab.
3. Scroll down to the Specimen Interpretation Factors card.
4. Click the Adequacy Code lookup field.
5. Click the Magnifying Glass icon or press Enter on your keyboard.
6. Select the updated code.
7. Click the Save button.

Changing the adequacy code will result in a reinterpretation being run and a new results mailer being generated.