



Food and Drug Branch – Device Recalls

CALIFORNIA DEVICE RECALL INFORMATION SHEET

Baxter Healthcare Corporation Recalls Abacus TPN Calculation Software for Potential Medication Error

Recall Date	Product Description	Recalling Firm	Recall Reason
6/22/2022	ABACUS TPN (Total Parenteral Nutrition) Calculation software Product codes: 8300-0167 (Abacus V3.1 CE), 8300-0168 to (Abacus V3.3 SE), 8300-3393 (Abacus V3.3 ME). Commonly used with ExactaMix automated compounding devices (ExactaMix 1200 and ExactaMix 2400).	Baxter Healthcare Corporation Deerfield, Illinois	There is a potential risk of medication error when using specific Abacus Order Entry & Calculation software.

Recall Class	Product Identification	Distribution	Affected Dates
I	Product Code 8300-3391, UDI (01)05413765577345(10)3.3.2.1, Release date 20-Sep-2019; to Product Code 8300-3393, UDI (01)05413765577369(10)3.3.2.1, Release date 20-Sep-2019; Product Codes: 8300-0167, V3.1 CE, Release date 20-May-2013, 8300-0168, V3.1 SE, Release date 20-May-2013 to 8300-0193, V3.2 ME, Release date 17-Dec-2015, No UDIs	Nationwide	June 2022 and prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

