



Food and Drug Branch – Drug Recalls

**CALIFORNIA DRUG RECALL INFORMATION SHEET**

**Siemens Medical Solutions Recalls Acuson Ultrasound Imaging System and Software for Calculation Error during 2D Manual Trace**

Recall Date	Product Description	Recalling Firm	Recall Reason
11/12/2021	<b>ACUSON Sequoia Ultrasound Imaging System with Software</b> Version VA25A/B	<b>Siemens Medical Solutions USA, Inc.</b> Issaquah, Washington	Due to a calculation error in the measurement when using 2D trace (manual trace) tool. The trace circumference value is overestimated and may potentially result in misdiagnosis of a patient's condition.

Recall Class	Product Identification	Distribution	Affected Dates
II	System Serial Numbers: 800101-802972	47 Units in California	November 2021 and prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

